

ACME UK INC LIMITED

Urine Collection Bag



Unique and easy to contemplate design offers provision of periodic monitoring of urine excretion in ill patients on an hourly basis.

Sample Collection Port for extracting urine samples by syringe

Optional fitting with inlet tube stopper and clip for emergency stopover

Optional Thread fitting with discharge tube stopper

Bottom drainage outlet design with standard connector

Optional transparent burette shaped volume meter provided which helps in accurate measurement of urine excretion and also makes provisions for

sample collection.

The parts are inter-fitted in such a way to ensure zero chances of cross contamination.

There is an additional provision of a 2000 ml urine bag opening directly into the burette meter, thus ensuring auto over flow containment and convenience to empty measured volume of urine into the bag.

A long kink resistant and flexible PVC Tube is attached to the measuring unit offering ease of use.

